

# Experience Lumity.

## 90 Premium

Comfort in Echo  
Speech in Car  
Speech in Loud Noise  
Music  
Comfort in Noise  
Speech in Noise  
Calm Situation

Streaming Speech  
Streaming Music

AutoSense OS™ 5.0

SmartSpeech Technology

Speech Enhancer  
SpeechSensor  
StereoZoom 2.0  
Dynamic Noise Cancellation  
Motion Sensor Hearing<sup>1</sup>  
ActiveVent™ Receiver<sup>2</sup>  
Roger™ Technology

Comfort in Echo (manual)  
Speech in Loud Noise (manual)

Speech in 360° (manual)  
EchoBlock  
SoundRelax

DuoPhone  
Real Ear Sound  
UltraZoom

WindBlock  
SoundRecover2  
WhistleBlock

Tinnitus Balance<sup>3</sup>  
T-coil<sup>4</sup>  
20 Channels

Enjoy easier hearing in a variety of situations<sup>1,2,3</sup> with Phonak SmartSpeech™ Technology in all levels of Lumity

## 70 Advanced

Music  
Comfort in Noise  
Speech in Noise  
Calm Situation  
Streaming Speech

Streaming Music

StereoZoom 2.0  
Dynamic Noise Cancellation  
Motion Sensor Hearing<sup>1</sup>  
ActiveVent Receiver<sup>2</sup>  
Roger Technology

Speech in Loud Noise (manual)  
Speech in 360° (manual)  
SoundRelax

DuoPhone  
Real Ear Sound  
UltraZoom

WindBlock  
SoundRecover2  
WhistleBlock

Tinnitus Balance<sup>3</sup>  
T-coil<sup>4</sup>  
20 Channels

## 50 Standard

Comfort in Noise  
Speech in Noise  
Calm Situation

Streaming Speech  
Streaming Music

Motion Sensor Hearing<sup>1</sup>  
ActiveVent Receiver<sup>2</sup>  
Roger Technology

DuoPhone  
Real Ear Sound  
UltraZoom

WindBlock  
SoundRecover2  
WhistleBlock

Tinnitus Balance<sup>3</sup>  
T-coil<sup>4</sup>  
16 Channels

## 30 Essential

Speech in Noise  
Calm Situation  
Streaming Speech  
Streaming Music

ActiveVent Receiver<sup>2</sup>  
Roger Technology  
UltraZoom

WindBlock  
SoundRecover2  
WhistleBlock

Tinnitus Balance<sup>3</sup>  
T-coil<sup>4</sup>  
12 Channels


All levels of Phonak Lumity include:



Universal connectivity



Compatibility with myPhonak and myPhonak Junior



Compatibility with Roger and Phonak wireless accessories

<sup>1</sup> Only available in rechargeable Lumity hearing aids

<sup>2</sup> Only available for Audéo L-R/RT/RL

<sup>3</sup> Not available in Naida L-UP or Sky Lumity hearing aids

<sup>4</sup> Only available in Audéo L-RT, Naida L-SP/UP and Sky L-M/SP/UP hearing aids

# Phonak Lumity

## Features and benefits

<b>Adaptive Phonak Digital 2.0</b>	Phonak proprietary fitting formula for reduced listening effort in noise for mild to severe hearing loss groups <sup>4***</sup>
<b>AutoSense OS 5.0</b>	Scans the environment and creates up to 200 unique setting combinations to accommodate your unique listening environment
<b>Phonak SmartSpeech Technology</b>	A collection of hearing aid features that help you participate in conversations and provides easier hearing in many situations <sup>1</sup>
<b>SpeechSensor</b>	Optimises your ability to hear speech from the sides and back in noise
<b>Speech Enhancer</b>	Easier to hear soft speech or speech from a distance in quiet situations <sup>1</sup>
<b>StereoZoom 2.0</b>	Automatic smooth activation of the hearing aid microphones to focus on a speaker in front of you in a challenging noisy environment
<b>Dynamic Noise Cancellation</b>	Automatically reduces background noise around you so that it is easier to hear speech in front of you <sup>1</sup>
<b>Motion Sensor Hearing*</b>	Detects whether you are standing or moving and adjusts the microphones to optimise your ability to hear in a noisy environment
<b>Roger Technology</b>	RogerDirect allows for direct streaming from Roger microphones, for better hearing in noise and over distance <sup>2</sup>
<b>ActiveVent Receiver</b>	Intelligently balances speech clarity and listening comfort. Provides superior sound quality for streaming media <sup>3**</sup>
<b>UltraZoom</b>	Automatically helps you focus on a speaker directly in front of you in a noisy environment
<b>SoundRecover2</b>	Makes high-pitched sounds easier to hear
<b>myPhonak</b>	The all-in-one app allows you to customise and adjust your hearing aids. You can also connect with your clinician remotely
<b>Tap Control*</b>	Activate voice assistant application, answer and end phone calls, and pause and resume streaming media
<b>Universal Connectivity</b>	Automatic streaming and hands-free calls to both ears from any Bluetooth® enabled device

\* Only available in rechargeable Lumity hearing aids  
 \*\* Compared to a standard acoustic coupling in quiet  
 \*\*\* Compared to APD

The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Sonova AG is under license.

- 1 Appleton, J. (2020) AutoSense OS 4.0 - significantly less listening effort and preferred for speech intelligibility. *Phonak Field Study News*, retrieved from [www.phonak.com/evidence](http://www.phonak.com/evidence).
- 2 Thibodeau L.M. (2020) Benefits in Speech Recognition in Noise with Remote Wireless Microphones in Group Settings. *Journal of the American Academy of Audiology*, 31(6), 404–411.
- 3 Latzel, M. & Hobi S. (2021) Receiver with mechanical vent provides benefit of open and closed acoustics for better speech understanding in noise and naturalness of own voice perception. *Phonak Field Study News*, retrieved from [www.phonak.com/evidence](http://www.phonak.com/evidence).
- 4 Wright, A. (2020). Adaptive Phonak Digital 2.0: Next-level fitting formula with adaptive compression for reduced listening effort. *Phonak Field Study News*, retrieved from [www.phonak.com/evidence](http://www.phonak.com/evidence).